

VIBRO-DEAERATOR for VACUUM PYCNOMETER

Code : A491



Used to vibrate the pycnometer to speed up the process and help air bubbles move faster resulting in shorter specimen preparation time

STANDARDS

ASTM D2041 • EN 12697-5 • EN 13108

ORDERING INFORMATION

Item	Code
VIBRO-DEAERATOR for VACUUM PYCNOMETER	A491X00XU
DEVICE for CLAMPING PYCNOMETER to SIEVE SHAKER	A491P001H
VACUUM PYCNOMETER - 10 lt	A490X00XH

VACUUM PYCNOMETER - 2 lt	A490X02XH
VACUUM PYCNOMETER - 5 lt	A490X05XH
SPRING CLAMP for A-490	A490P001H

ORDERING INFORMATION