

ERLENMEYER

Code : W002



- An Erlenmeyer flask in engineering testing laboratories is used for mixing, heating, and storing liquids or samples. Its conical shape allows for easy swirling without spilling, and it is commonly used in titrations, chemical reactions, and sample preparation, offering both stability and ease of use.

ORDERING INFORMATION

Item Name	Item Code
ERLENMEYER - 1000 ml	W002G01KH
ERLENMEYER - 25 ml	W002G025H
ERLENMEYER - 2000 ml	W002G02KH
ERLENMEYER - 3000 ml	W002G03KH

ERLENMEYER - 50 ml W002G050H

ERLENMEYER - 5000 ml W002G05KH

ERLENMEYER - 100 ml W002G100H

ERLENMEYER - 150 ml W002G150H

ERLENMEYER - 250 ml W002G250H

ERLENMEYER - 300 ml W002G300H

ERLENMEYER - 500 ml W002G500H

OTHER PHOTOS

