

GLASS FUNNEL

Code : W020



- A glass funnel in engineering testing laboratories is used for transferring liquids or powders into containers with narrow openings. It helps ensure a smooth, controlled flow, preventing spills and contamination, and is commonly used in filtration, sample preparation, and liquid handling processes.

ORDERING INFORMATION

Item Name	Item Code
GLASS FUNNEL - Ø 35 mm	W020G006H
GLASS FUNNEL - Ø 55 mm	W020G008H
GLASS FUNNEL - Ø 80 mm	W020G010H
GLASS FUNNEL - Ø 120 mm	W020G013H
GLASS FUNNEL - Ø 150 mm	W020G016H

GLASS FUNNEL - Ø 200 mm

W020G026H