

FILTERING FLASK

Code : W021



- A filtering flask in engineering testing laboratories is used to hold liquids during filtration processes. It is typically used with a vacuum pump to create suction, which speeds up the filtration of liquids through a filter paper, ensuring the separation of solids from liquids in various chemical or material testing procedures.

ORDERING INFORMATION

Item Name	Item Code
FILTERING FLASK - 1000 ml	W021G01KH
FILTERING FLASK - 2000 ml	W021G02KH
FILTERING FLASK - 100 ml	W021G100H

FILTERING FLASK - 250 ml

W021G250H

FILTERING FLASK - 500 ml

W021G500H

OTHER PHOTOS

