LAB TOOLS > GLASSWARE

GLASS BEAKER

Code : W001



- A glass beaker in engineering testing laboratories is used for holding, mixing, heating, and measuring liquids or powders. Its transparent nature allows for easy observation, and it is commonly used in chemical reactions, sample preparation, and general laboratory tasks.

ORDERING INFORMATION

Item Name	Item Code
GLASS BEAKER - 5 ml	W001G005H
GLASS BEAKER - 10 ml	W001G010H
GLASS BEAKER - 1000 ml	W001G01KH
GLASS BEAKER - 25 ml	W001G025H

GLASS BEAKER - 2000 ml	W001G02KH
GLASS BEAKER - 3000 ml	W001G03KH
GLASS BEAKER - 50 ml	W001G050H
GLASS BEAKER - 5000 ml	W001G05KH
GLASS BEAKER - 100 ml	W001G100H
GLASS BEAKER - 150 ml	W001G150H
GLASS BEAKER - 250 ml	W001G250H
GLASS BEAKER - 400 ml	W001G400H
GLASS BEAKER - 600 ml	W001G600H
GLASS BEAKER - 800 ml	W001G800H

OTHER PHOTOS

