LAB TOOLS > GLASSWARE

# FILTERING FLASK

Code : W021



- A filtering flask in engineering testing laboratories is used to hold liquids during filtration processes. It is typically used with a vacuum pump to create suction, which speeds up the filtration of liquids through a filter paper, ensuring the separation of solids from liquids in various chemical or material testing procedures.

## ORDERING INFORMATION

| Item Name | Item Code |
|---|---|
| FILTERING FLASK - 1000 ml | W021G01KH |
| FILTERING FLASK - 2000 ml | W021G02KH |
| FILTERING FLASK - 100 ml | W021G100H |

| | |
|---|---|
| FILTERING FLASK - 250 ml | W021G250H |
| FILTERING FLASK - 500 ml | W021G500H |

## OTHER PHOTOS

