LAB TOOLS > GLASSWARE

PORCELAIN MORTAR and PESTLE

Code : W026



- A porcelain mortar and pestle in engineering testing laboratories are used for grinding, crushing, and mixing solids into fine powders. The non-porous nature of porcelain ensures that no material is absorbed, making it ideal for preparing samples for chemical or material analysis without contamination.

ORDERING INFORMATION

Item Name	Item Code
PORCELAIN MORTAR and PESTLE - Ø 80 mm	W026C080H
PORCELAIN MORTAR and PESTLE - Ø 100 mm	W026C100H
PORCELAIN MORTAR and PESTLE - Ø 130 mm	W026C130H
PORCELAIN MORTAR and PESTLE - Ø 160 mm	W026C160H

PORCELAIN MORTAR and PESTLE - Ø 190 mm W026C190H

OTHER PHOTOS

