PLASTIC WASHING BOTTLE

Code : W031



- A washing bottle in engineering testing laboratories is used for rinsing or cleaning laboratory equipment, such as beakers, flasks, or test tubes. It typically contains distilled or deionized water and allows for controlled and directed rinsing to remove residues or contaminants from samples or instruments.

ORDERING INFORMATION

Item Name	Item Code
WASHING BOTTLE - 1000 ml	W031P01KH
WASHING BOTTLE - 100 ml	W031P100H
WASHING BOTTLE - 250 ml	W031P250H
WASHING BOTTLE - 500 ml	W031P500H